

Acrylic RPD Options

						
	Acetal Tooth Shade Framework	Acrylic	Clear Framework	CR / CO Cast Framework	Titanium Framework	Valplast
Comfort	Excellent	Moderate	Excellent (Best Choice)	Great	Great	Good
Aesthetics	Great	Moderate	Excellent	Good	Good	Excellent
Stability/Fit	Excellent	Good	Excellent (Best Choice)	Great	Great	Good
Weight	Light	Heavy	Light	Heavy	Moderate (3x lighter than CR/CO)	Moderate
Durability	Great	Moderate	Great	Excellent	Excellent	Great
Reliability/Repairability	Moderate	Easy	Moderate	Moderate	Moderate	Difficult
Dentition Friendly	Excellent	Moderate	Excellent	Moderate	Good	Good
Clasp Type	Tooth Shade A1 – D2	Wire Metal	Clear	Metal	Titanium	Valplast Gingival
Manufacturing Method	 Milled	 Analog	 Milled	 Cast	 Milled	 Analog